

FUJIFILM

Wako

Code No. 011-21891 (100 μ L)

**Anti Glutathione S-transferase, Monoclonal Antibody,
Peroxidase Conjugated**
**抗グルタチオン S - トランスフェラーゼ,
モノクローナル抗体, ペルオキシダーゼ結合**

Glutathione S-transferase (GST) has been used as the tag for detection and/or purification of recombinant protein.

Monoclonal Anti GST antibody was produced by hybridoma clone 5A7. This clone was established by fusing mouse myeloma cells and spleen cells of mice immunized with GST protein. This antibody was purified from mouse ascites and was conjugated with horseradish peroxidase. This product is for laboratory use only ; use in any such application is the responsibility of the user.

【Concentration】

Indicated on the label.

【Formulation】

1 \times PBS aqueous solution with 50% glycerol.

【Clone No.】

5A7

【Subclass】

IgG_{2b}

【Specificity】

Specific for Glutathione S-transferase.

【Working Dilution】

Western blot 1 : 250 – 1 : 1,000

【Storage】

Keep at -20°C in the dark.

Avoid repeated freeze and thaw.

【Package】

Code No.	Packaging
011-21891	100 μ L

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2304KA2