

FUJIFILM

Wako

<for Research Use Only> Code No. 015-24231 (10 mg)
011-24233 (50 mg)

Avidin, Neutralized, from Egg White
中性アビジン, 卵白由来

Neutralized avidin is a chemically and enzymatically modified avidin. In the result of these modifications, Neutralized avidin has been lower in non-specific binding than normal avidin from egg white. Furthermore, the isoelectric point (pI) has been shifted to the neutral pH.

Source : Egg white

Appearance : Lyophilized

Molecular weight : approximately 14.7 kDa as monomer.

Biotin binding capacity : Indicated on the label.

Unit Definition : One unit is the binding quantity of 1 μ g of biotin.

[Storage]

Store at -20°C in the dark.

[Package]

Code No.	Packaging
015-24231	10 mg
011-24233	50 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggenstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA2