

FUJIFILM

Wako

Code No. 012-22801 (50 μ L)

**Anti Human Diacylglycerol Kinase δ ,
Rabbit (for Immunoprecipitation)**
抗ヒトジアシルグリセロールキナーゼ δ ,
ウサギ(免疫沈降用)

Rabbit anti human DGK δ antibody is raised against synthetic peptide corresponding to C-terminals of δ 1 or δ 2 of human DGK δ .

This product is for laboratory use only; use in any such application is the responsibility of the user.

Preparation : Purified by protein A affinity chromatography from rabbit antisera and prepared in PBS solution with 50% glycerol containing no preservatives and stabilizers.

Specificity : Specific to human DGK δ

Working Quantity : Immunoprecipitation 5~10 μ g/test

Storage : Keep at -20°C

Package : 50 μ L

Concentration : Indicated on the label.

Reference :

1) Imai, S. *et al.* : *Biochim. Biophys. Acta*, **1791**, 246 (2009).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA2