

FUJIFILM

Wako

Code No. 012-27541 (100 μ L)

Anti Sox2, Rabbit

Sox2 is a member of the SRY (sex determining region-Y)-related HMG-box family of transcriptional factors, and is expressed in inner cell mass, epiblast and extraembryonic mesoderm. Sox2 is required to maintain the undifferentiated state of pluripotent stem cell. Sox2 is reported as one of the factors to generate induced pluripotent stem cells.

This product is for laboratory research use only ; use in any such application is the responsibility of the user.

[Formulation] :

Phosphate buffer containing 50w/v% Glycerol

[Source] :

Purified from rabbit serum using affinity column

[Application] :

Western Blotting 1 : 1,000

Immunocytochemistry 1 : 200~500

※Please optimize the appropriate concentration for your analysis conditions and samples.

[Package] :

100 μ L

[Storage] :

Store at -20°C .

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA2