

**FUJIFILM**

**Wako**

Code No. 013-21851 (200  $\mu$ g)  
017-21854 ( 1 mg)  
013-21856 ( 5 mg)

**Anti Glutathione S-transferase, Monoclonal Antibody**  
**抗グルタチオン S - トランスフェラーゼ ,**  
**モノクローナル抗体**

Glutathione S-transferase (GST) has been used as the tag for detection and/or purification of recombinant protein.

Monoclonal Anti GST antibody was produced by hybridoma clone 5A7. This clone was established by fusing mouse myeloma cells and lymphocytes of mice immunized with GST protein. This antibody was purified from mouse ascites.

This product is for laboratory use only ; use in any such application is the responsibility of the user.

**【Concentration】**

Indicated on the label.

**【Formulation】**

1  $\times$  PBS aqueous solution with 50% glycerol.

**【Clone No.】**

5A7

**【Subclass】**

IgG<sub>2b</sub>

**【Specificity】**

Specific for Glutathione S-transferase.

**【Working Concentration】**

Western blot 0.2-2  $\mu$ g/mL

Immunoprecipitation 1-2  $\mu$ g/20  $\mu$ L Protein A or Protein G beads slurry

**【Storage】**

Keep at -20°C in the dark.

Avoid repeated freeze and thaw.

**【Package】**

Code No.	Packaging
013-21851	200 $\mu$ g
017-21854	1 mg
013-21856	5 mg

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-31110  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

2304KA2