

FUJIFILM**Wako**Code No. 018-26921 (10 μ L)
014-26923 (50 μ L)**Anti Amyloid β 40, Monoclonal Antibody (BA27)**
抗アミロイド β 40, モノクローナル抗体 (BA27)**Background**

Amyloid β (A β) is a peptide which consists of about 40 amino acids. A β is considered the cause of Alzheimer's disease, because it is a major component of senile plaques formed in the brains of Alzheimer's disease patients. A β has several types differing in length of the amino acid, The main types are A β 40 (40 amino acid), A β 42 (42 amino acid), and A β 43 (43 amino acid). This item is a mouse monoclonal antibody that specifically reacts with A β 40¹⁾. The antibody was developed by Takeda Pharmaceutical Company Limited.

Description

[Reactivity]	Reacts with Amyloid β 40
[Antigen]	1-40 a.a of human Amyloid β
[Species cross reactivity]	Human, Mouse and Rat
[Clone No.]	BA27
[Host]	Mouse
[Isotype]	IgG2a $\cdot$ κ
[Purification]	Protein A-affinity purified
[Concentration]	Indicated on the label
[Formulation]	PBS with 0.05% Sodium Azide

Applications

Western Blot ²⁾³⁾⁶⁾	1 : 100-1,000
Immunohistochemistry ⁴⁾⁵⁾	1 : 500-5,000
ELISA ¹⁾²⁾³⁾⁶⁾	1 : 100-1,000

Optimal concentration should be determined by each laboratory for each application.

Storage

Store at -20 °C.
Avoid repeated freeze and thaw.

Package

Code No.	Packaging
018-26921	10 μ L
014-26923	50 μ L

References

- 1) Asami-Odaka, A., *et. al.* : *Biochemistry*, **34**, 32, 10272 (1995).
- 2) Tomita, T., *et. al.* : *J. Neurosci.* **19**, 24, 10627 (1999).
- 3) Sudoh, S., *et. al.* : *J. Neurochem.* **71**, 4, 1535 (1998).
- 4) Iwatsubo, T., *et. al.* : *Am J. Pathol.*, **149**, 6, 1823 (1996).
- 5) Nunomura, A., *et al.* : *J. Neuropathol. Exp. Neurol.*, **59**, 11, 1011 (2000).
- 6) Iwata, H., *et. al.* : *J. Biol. Chem.*, **276**, 24, 21678 (2001).

Acknowledgement

The Purchaser shall use the Product at its research facility solely for its internal research. A separate license by Takeda Pharmaceutical Company Limited (1-1, Doshomachi 4-chome, Chuo-ku, Osaka 540-8645, Japan) is required for Commercial Use of the Product. Information about a license for the Commercial Use of the Product may be obtained from :

FUJIFILM Wako Pure Chemical Corporation. Laboratory Chemicals Division, 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>