

FUJIFILM

Wako

Code No. 015-22391 (200 μ L)
019-22394 (1 mL)

**Anti DYKDDDDK tag, Monoclonal Antibody,
Peroxidase Conjugated**

**抗DYKDDDDKタグ, モノクローナル抗体
ペルオキシダーゼ結合**

The DYKDDDDK was consisted of 8 amino acids and has been used as the tag for detection and/or purification of recombinant protein.

Monoclonal anti DYKDDDDK tag antibody was established by fusing mouse myeloma cells and spleen cells of mice immunized with synthetic DYKDDDDK peptide which was coupled with keyhole limpet hemocyanin (KLH) to mouse. This antibody was purified from mouse ascites.

This antibody was conjugated with horse radish peroxidase.

[Concentration]

Describing at label.

[Formulation]

1 $\times$ PBS aqueous solution with 50% glycerol.

[Clone No.]

1E6

[Subclass]

IgG_{2b}

[Specificity]

Specific for the DYKDDDDK peptide.

[Working Dilution]

Western blot 1 : 1,000 – 1 : 20,000

[Storage]

Keep at -20°C .

Do not freeze and thaw cycle.

[Package]

Code No.	Packaging
015-22391	200 μ L
019-22394	1 mL

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwkw.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA2