

FUJIFILM**Wako**Code No. 011-23591 (200 μ g)
015-23594 (1 mg)**Anti V5 tag, Monoclonal Antibody (6F5)**
抗V5タグ, モノクローナル抗体(6F5)

The GKPIPNNLLGLDST peptide derived from Simian Virus 5 was consisted of 14 amino acids and has been used as the tag for detection and/or purification of recombinant protein. Monoclonal anti V5 antibody was produced by hybridoma clone 6F5. This clone was established by immunization of synthetic GKPIPNNLLGLDST peptide which was coupled with keyhole limpet hemocyanin (KLH) to mouse. This antibody was purified by affinity chromatography on Protein A from mouse ascites. This product is for laboratory use only ; use in any such application is the responsibility of the user.

[Concentration]

Indicated on the label.

[Formulation]

1 × PBS aqueous solution with 50% glycerol.

[Clone No.]

6F5

[Subclass]IgG_{2b}**[Specificity]**

Specific for the GKPIPNNLLGLDST peptide.

[Working Dilution]Western blot 1 : 1,000 – 1 : 10,000
Immunoprecipitation 2 – 10 μ g/20 μ L 10% Protein A or Protein G beads slurry

Please optimize the most appropriate concentration for your analysis conditions and samples.

[Storage]Keep at – 20 °C in the dark.
Do not freeze and thaw cycle.**[Package]**

Code No.	Packaging
011-23591	200 μ g
015-23594	1 mg

FUJIFILM Wako Pure Chemical Corporation1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>**FUJIFILM Wako Chemicals U.S.A. Corporation**1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>**FUJIFILM Wako Chemicals Europe GmbH**Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2