

FUJIFILM**Wako**

Code No. 015-25691 (100 µg)

Anti Mutated IDH1/2, Monoclonal Antibody (MsMab-1)**抗変異型 IDH1/2, モノクローナル抗体 (MsMab-1)****Background**

Mutations of isocitrate dehydrogenase 1 (IDH1) and isocitrate dehydrogenase 2 (IDH2) were found to create the capability of the enzyme to change α -ketoglutarate to an oncometabolite, R(-)-2-hydroxyglutarate (2-HG), although IDH1 and IDH2 catalyze the oxidative carboxylation of isocitrate to α -ketoglutarate in cytosol and mitochondria, respectively. IDH1/2 mutations occur in gliomas, cartilaginous tumors, osteosarcomas, and in acute myeloid leukemias. In astrocytomas, oligodendrogliomas, oligoastrocytomas, and secondary glioblastomas, IDH1/2 mutations have been identified as early and frequent genetic alterations (50-93%). They might be the initiating event in these glioma subtypes. In contrast, primary glioblastomas rarely contain IDH1/2 mutations (<5%). The IDH1 mutations are remarkably specific to a single codon in the conserved and functionally important arginine 132 residue (R132). In contrast, the IDH2 mutations are specific to a single codon in arginine 172 residue (R172) in gliomas. The vast majority of changes are heterozygous. In gliomas, IDH1 mutations were reported as IDH1-R132H (664/716 : 92.7%), IDH1-R132C (29/716 : 4.2%), IDH1-R132S (11/716 : 1.5%), IDH1-R132G (10/716 : 1.4%), and IDH1-R132L (2/716 : 0.2%); IDH2 mutations were reported as IDH2-R172K (20/31 : 64.5%), IDH2-R172M (6/31 : 19.3%), and IDH2-R172W (5/31 : 16.2%).

This Anti Mutated IDH1/2, Monoclonal Antibody is available to detect IDH1-R132H/R132S/R132G, IDH2-R172M/R172S/R172G.

This product is for laboratory use only ; use in any such application is the responsibility of the user.

Description

[Reactivity]	IDH1-R132H/R132S/R132G, IDH2-R172M/R172S/R172G
[Species cross reactivity]	-
[Clone No.]	MsMab-1
[Host]	Mouse
[Isotype]	IgG _{2a} , kappa
[Source]	Purified from hybridoma supernatant using protein G
[Concentration]	Indicated on the label
[Formulation]	PBS (pH 7.4), containing 0.1% sodium azide

Applications

ELISA	1 µg/mL
Western blot	1 µg/mL
Immunohistochemistry	5 µg/mL

Storage

Store at -20 °C.
Avoid repeated freeze and thaw.

- 1/2 -

Package

100 µg

References

- 1) Kaneko, M. K., Ogasawara, S. and Kato, Y. : *TJEM*, **230**, 103 (2013).
- 2) Liu, X., Kato, Y., Kaneko, M. K., Sugawara, M., Ogasawara, S., Tsujimoto, Y., Naganuma, Y., Yamakawa, M., Tsuchiya, T. and Takagi, M. : *Cancer Med.*, **2**(6), 803 (2013).
- 3) Ogasawara, S., Kaneko, M. K., Tsujimoto, Y., Liu, X. and Kato, Y. : *Monoclon. Antib. Immunodiagn. Immunother.*, **32**(6), 377 (2013).
- 4) Moriya, K., Kaneko, M. K., Liu, X., Hosaka, M., Fujishima, F., Watanabe, M., Ogasawara, S., Sakuma, J., Sasahara, Y., Kure, S. and Kato, Y. : *Cancer Sci.*, in press

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2

- 2/2 -