

FUJIFILM

Wako

Code No. 015-27531 (100 μ L)

Anti Oct3/4, Rabbit

Oct3/4 is a member of the POU family of transcription factors, and is expressed in embryonic stem cells and embryonic carcinoma cells. Oct3/4 regulates the expression of Tc11 (T cell leukemia/lymphoma 1) and Nanog, and contributes to cell proliferation and stabilization of cell pluripotency. Oct3/4 is reported as one of the factors to generate induced pluripotent stem cells.

This product is for laboratory research use only ; use in any such application is the responsibility of the user.

[Formulation] :

Phosphate buffer containing 50w/v% Glycerol

[Source] :

Purified from rabbit serum using affinity column

[Application] :

Western Blotting 1 : 1,000

Immunocytochemistry 1 : 500

※ Please optimize the appropriate concentration for your analysis conditions and samples.

[Package] :

100 μ L

[Storage] :

Store at -20°C .

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA2