

FUJIFILM

Wako

Code No. 016-23041 (500 μ g)

**Apolipoprotein A-I Solution,
from Human Plasma**
アポリポプロテインA-I 溶液, ヒト血漿由来

Source : Human plasma

Formulation : Liquid. 10 mmol/L Ammonium Hydrogen Carbonate, pH 7.4

Molecular weight : approx. 28,000

Protein concentration : Indicated on the label

[Expiration date] Indicated on the label.

[Storage] Store at -80°C .

[Package] 500 μ g

[Caution] Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.
This product is purified from human plasma which is tested and found negative for HBsAg, anti-HCV, anti-HBc and negative for anti-HIV 1 & 2. However it should be handled as potentially infectious material.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA2