

FUJIFILM

Wako

Code No. 018-22141 (100 μ L)

Anti SQSTM1/A170/p62, Rabbit 抗SQSTM1/A170/p62, ウサギ

Anti SQSTM1/A170/p62 antibody is raised in rabbits using recombinant murine SQSTM1/A170/p62 PEST domain (T7 tag + amino acid 254-333 + His tag) as the immunogen.

Preparation : The antiserum is diluted with an equal amount of PBS and absorbed with SQSTM1/A170/p62 knockout mouse tissue proteins.

Formulation : Diluted antiserum containing 0.1 % Sodium azide

Specificity : Specific for mouse and rat SQSTM1/A170/p62. Slightly reactive with human SQSTM1/A170/p62.

Working dilution :

Western blot	1 : 200 ~ 1 : 1,000
Immunohistochemistry	1 : 400 ~ 1 : 1,000
Immunofluorescence	1 : 400 ~ 1 : 1,000

[Storage] Store at -20°C .

[Package] 100 μ L

[Reference]

- 1) Ishii, T., Itoh, K., Takahashi, S., Sato, H., Yanagawa, T., Katoh, Y., Bannai, S. and Yamamoto, M. : *J. Biol. Chem.*, **275**, 16023 (2000).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2304KA2