

FUJIFILM**Wako**Code No. 018-24101 (100 μ g)**Anti Human Podoplanin, Monoclonal Antibody**
抗ヒトポドプランニン, モノクローナル抗体**Background**

Podoplanin, also known as Aggrus, gp44, or T1 α , is a type I transmembrane protein, which possesses platelet-aggregating activity and metastatic activity. In addition to this, Podoplanin is known as a lymphatic endothelium-specific marker, and is over-expressed in several cancers. This Anti Human Podoplanin, Monoclonal Antibody is available to detect human podoplanin. This product is for laboratory use only ; use in any such application is the responsibility of the user.

Description

[Reactivity]	Reacts with human podoplanin
[Species cross reactivity]	–
[Clone No.]	NZ-1.2
[Host]	Rat
[Isotype]	IgG2a
[Source]	Purified from rat ascites
[Concentration]	Indicated on the label
[Formulation]	PBS (pH 7.2) containing 0.05% sodium azide

Applications

ELISA	0.5-1	μ g/mL
Western blot	0.1-1	μ g/mL
Flowcytometry	1-5	μ g/mL
Immunoprecipitation	1-5	μ g/mL
Immunohistochemistry	1-5	μ g/mL

Storage

Store at -20°C .
Avoid repeated freeze and thaw.

Package100 μ g**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggenstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2304KA2