

FUJIFILM

Wako

Code No. 018-27521 (100 μ L)

Anti Nanog, Rabbit

Nanog is a transcription factor and is expressed early embryos and pluripotent stem cells including mouse and human embryonic stem cells. Nanog is required for maintenance of pluripotency.

This product is for laboratory research use only ; use in any such application is the responsibility of the user.

[Formulation] :

Phosphate buffer containing 50w/v% Glycerol

[Source] :

Purified from rabbit serum using affinity column

[Application] :

Western Blotting 1 : 500

Immunocytochemistry 1 : 200

※Please optimize the appropriate concentration for your analysis conditions and samples.

[Package] :

100 μ L

[Storage] :

Store at -20°C .

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA2