

FUJIFILM

Wako

Code No. 034-26241 (20 μ g)
030-26243 (1 mg)

CTGF, Human, recombinant
結合組織成長因子, ヒト, 組換え体
(rhCTGF)

Source : *E. coli* expressed human CTGF

Formulation : Lyophilized from a 0.2 μ m filtered solution containing 0.1% Trifluoroacetic Acid.

Molecular weight : 11,200

Biological activity : ED₅₀ is less than 2000 ng/mL determined by the HUVEC cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other aqueous buffers.

Endotoxin : Less than 0.1 EU/ μ g

[Expiration date]

Indicated on the label.

[Storage]

Store at -20°C.

[Package]

Code No.	Packaging
034-26241	20 μ g
030-26243	1 mg

[Caution]

Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors

