

FUJIFILM

Wako

Code No. 031-25911 (100 mg)
037-25913 (250 mg)
035-25914 (1 g)

Concanavalin A コンカナバリン A

Concanavalin A is a plant-derived lectin. It is reported that agglutinates blood cells and precipitates polysaccharides and glycoproteins. It also acts as a mitogen for lymphocytes.

Source : Sword Beans

Appearance : White ~ brown, Crystals ~ crystalline powder

Reconstitution (example) : Reconstitute in 10 mmol/L Phosphate Buffer (pH 7.0), 0.15 mol/L NaCl to a concentration of 0.2 mg/mL. *Do not vortex.*

Carbohydrate specificity : α -Mannose, α -Glucose

【Expiration data】

Indicated on the label.

【Storage】

Store in the dark at -20°C.

【Package】

Code No.	Package
031-25911	100 mg
037-25913	250 mg
035-25914	1 g

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://flwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA1