

FUJIFILM**Wako**Code No. 038-21543 (5 mg)
036-21544 (25 mg)

COA-Cl

COA-Cl is found to be an angiogenesis promoter.
This promoting action is expressed by the phosphorylation of
one MAPK-signaling factor.


Appearance : White ~ pale brown, crystalline powder or lump

Solubility : Soluble in water, PBS and normal saline

Molecular Formula : C₁₁H₁₄ClN₅O₂

Molecular Weight : 283.71

Structure :



[Expiration Date] Indicated on the label.

[Storage] Store at room temperature.
After reconstitution, store at 2-8°C

[Packages]

Code No.	Packaging
038-21543	5 mg
036-21544	25 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://flwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>