

FUJIFILM

Wako

Code No. 047-34581 (1 mg)
043-34583 (5 mg)
041-34584 (25 mg)

DYKDDDDK Peptide **DYKDDDDKペプチド**

The DYKDDDDK Peptide is used for the competitive elution of DYKDDDDK tag fusion proteins from the immune complex after immunoprecipitation by using anti DYKDDDDK tag antibody. This product is for laboratory use only ; use in any such application is the responsibility of the user.

【Formulation】

Trifluoroacetic acid salt.

【Amino acid sequence】

Asp-Tyr-Lys-Asp-Asp-Asp-Lys

【Preparation of stock solution】

Dissolve in 1×TBS (10 mmol/L Tris-HCl, pH 7.4, 150 mmol/L Sodium Chloride) to a final concentration of 5 mg/mL. Aliquot and store at -20°C. Avoid repeated freeze and thaw cycle.

【Working concentration】

Standard : 100 µg/mL in 1×TBS

A working concentration depends on the experiment.

【Storage】

Store at -20°C and protect from light.

【Package】

Code No.	Packaging
047-34581	1 mg
043-34583	5 mg
041-34584	25 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://flwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA3