

FUJIFILM

Wako

< for Research Use Only >

Code No. 042-33431 (0.5 mg)

**Discoidin II Lectin, recombinant,
Solution**

Discoidin II レクチン, 組換え体, 溶液

Discoidin II lectin is isolated from *Dictyostelium discoideum*. It is known that binds LacNAc and galactose.

The slime mold *Dictyostelium discoideum* produces both Discoidin I and Discoidin II. The localization of Discoidin I is prominent in aggregating cells, and Discoidin II is in fruiting body formation and is localized in the prespore vesicles.

Source : *E. coli* expressed *Dictyostelium discoideum* lectin (Discoidin II).

Appearance : PBS (pH 7.4).

Protein concentration : Indicated on the label.

Glycan specificity : Asialoglycans. LacNAc. Gal

Manufacture : Affinity purification by lactose agarose.

[Expiration date] Indicated on the label.

[Storage] Store at -20°C . After thawing, keep at $2-10^{\circ}\text{C}$. Avoid the freeze-thaw cycles.

[Package] 0.5 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2