

FUJIFILM

Wako

Code No. 058-09521 (100 μ g)
054-09523 (500 μ g)

**Epidermal Growth Factor,
Human, recombinant, Animal-derived-free**
上皮細胞成長因子,
ヒト, 組換え体, 動物由来物フリー
(rhEGF, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human epidermal growth factor

Formulation : Lyophilized from a 0.2 μ m filtered solution containing 0.1% Trifluoroacetic Acid.

Molecular weight : 6,200 (calculated)

Biological activity : ED₅₀ is less than 0.1 ng/mL determined by the BALB/c 3T3 cell proliferation.

Reconstitution : Reconstitute in sterile water to a concentration of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other aqueous buffers.

Endotoxin : Less than 0.1 EU/ μ g.

【Expiration date】 Indicated on the label.

【Storage】 Store at -20°C .

【Package】

Code No.	Packaging
058-09521	100 μ g
054-09523	500 μ g

【Caution】 Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2309KA1