

FUJIFILM

Wako

Code No. 064-07081 (10 μ g)
060-07083 (1 mg)

**Flt3 Ligand, Human, recombinant,
Animal-derived-free**
Flt3 リガンド, ヒト, 組換え体, 動物由来物フリー
(rhFlt3-L, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human Flt3 ligand

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 10 mmol/L Sodium Phosphate, 50
mmol/L Sodium Chloride, pH 7.5.

Molecular weight : 17,700 (calculated)

Biological activity : ED₅₀ is less than 5 ng/mL determined by
the OCI-AML5 cell proliferation.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*

This solution can then be diluted into other
aqueous buffers.

Endotoxin : Less than 0.1 EU/ μ g

[Expiration date] Indicated on the label.

[Storage] Store at -20°C.

[Package]

Code No.	Packaging
064-07081	10 μ g
060-07083	1 mg

[Caution] Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://flwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2309KA1