

FUJIFILM**Wako**


Code No. 061-06273 (5 mg)

Fluorescein Di- β -D-galactopyranoside (FDG)

The Fluorescein Di- β -D-galactopyranoside (FDG) is a fluorescence probe that reacts to the enzyme β -galactosidase in cells. This product is for laboratory use only ; use in any such application is the responsibility of users.

Appearance : White~Nearly white, Crystals~Powder

Structure :



Molecular formula : $C_{32}H_{32}O_{15}$

Molecular weight : 656.59

Solubility : Methanol

Ex/Em : 490 ± 10 nm/ 515 ± 10 nm
(in Methanol)

Storage : Store away from sunlight at -20°C ,
under inert gas.

Package :

Code No.	Packaging
061-06273	5 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggenstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2