

FUJIFILM

Wako

< for Research Use Only >

Code No. 062-06281 (1 mg)

F17AG Lectin, recombinant, Solution

F17AGレクチン，組換え体，溶液

F17AG lectin is isolated from *E. coli* fimbriae. It is known that specifically bind to GlcNAc.

Source : *E. coli* expressed *E. coli* lectin (F17AG).

Appearance : 0.1 × PBS (pH 7.4).

Protein concentration : Indicated on the label.

Carbohydrate specificity : GlcNAc

Manufacture : Affinity purification by *N*-acetyl glucosamine agarose.

[Expiration date] Indicated on the label.

[Storage] Store at −20 °C. After thawing, keep at 2-10 °C. Avoid the freeze-thaw cycles.

[Package]

Code No.	Packaging
062-06281	1 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2