

FUJIFILM

Wako

Code No. 066-06681 (1 g)
062-06683 (5 g)

Fibrinogen, from Bovine Plasma フィブリノーゲン, ウシ血漿由来

Source : Bovine Plasma

Appearance : White to slightly brown, crystalline powder to powder

Formulation : Lyophilized from 20 mmol/L Citric Acid, 100 mmol/L NaCl (pH 6.8)

Clottable Protein : Indicated on the label

[Caution]

- 1) Solution may contain a small amount of insoluble aggregated material. Do not stir intensely or sonicate. Stirring or sonication may cause an increase or may generate additional aggregate.
- 2) To remove insoluble aggregated material, filter the fibrinogen solution.
- 3) Moisture is easily absorbed, so after using, please cap or seal immediately.

[Storage] Store at 2~10 °C with inert gas.

[Expiration date] Indicated on the label

[Package]

Code No.	Packaging
066-06681	1 g
062-06683	5 g

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA3