

FUJIFILM

Wako

Code No. 069-07075 (10 μ g)
067-07071 (50 μ g)
063-07073 (100 μ g)
061-07074 (1 mg)

**Fibroblast Growth Factor(basic), Human,
recombinant, Animal-derived-free**
線維芽細胞成長因子(塩基性), ヒト, 組換え体,
動物由来物フリー
(rh basic FGF/FGF2, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human fibroblast growth factor (basic).

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 10 mmol/L Sodium Phosphate, 75
mmol/L Sodium Chloride, pH 7.5.

Molecular weight : 17,300

Biological activity : ED₅₀ is below 0.5 ng/mL determined by
NR6R-3T3 cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 0.1 EU/ μ g

【Expiration date】

Indicated on the label.

【Storage】

Store at -20°C.

【Package】

Code No.	Packaging
069-07075	10 μ g
067-07071	50 μ g
063-07073	100 μ g
061-07074	1 mg

【Caution】

Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors

