

FUJIFILM

Wako

Code No. 062-06661 (50 μ L)
068-06663 (50 μ L $\times$ 4)

bFGF Solution, MF **bFGF 溶液, MF**

Source : *E. coli* expressed human fibroblast growth factor (basic)

Formulation : Sterilized solution containing Sucrose and
Disodium Edetate Hydrate

Specific Activity : $\geq 500,000$ IU/mg

Endotoxin : ≤ 10 EU/mg

Storage : Store at -20°C

Reference : Akutsu, H., Machida, M., Kanzaki, S., Sugawara,
T., Ohkura, T., Nakamura, N., Yamazaki, M.,
Yamazaki-Inoue, M., Miura, T., Vemuri, CM.,
Rao, SM., Miyado, K. and Umezawa, A. :
Regenerative Therapy, **1**, 18 (2015).

Package :

Code No.	Packaging
062-06661	50 μ L
068-06663	50 μ L $\times$ 4

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2