

**FUJIFILM**

**Wako**

< for Research Use Only >

Code No. 070-06381 (1 mg)

**GRFT Lectin, recombinant, Solution**  
**GRFTレクチン, 組換え体, 溶液**

GRFT is isolated from *Griffithia* sp., the red algae. It is known that bind potently to high-mannose structure and be reported the inhibition of virus infection.

**Source** : *E. coli* expressed *Griffithia* sp. lectin (GRFT).

**Appearance** : PBS (pH 7.4).

**Protein concentration** : Indicated on the label.

**Carbohydrate specificity** : High-mannose

**Manufacture** : Affinity purification by mannose agarose.

[**Expiration date**] Indicated on the label.

[**Storage**] Store at  $-20^{\circ}\text{C}$ . After thawing, keep at  $2-10^{\circ}\text{C}$ .  
Avoid the freeze-thaw cycles.

[**Package**]

| Code No.  | Packaging |
|-----------|-----------|
| 070-06381 | 1 mg      |

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-311-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

2305KA2