

FUJIFILM

Wako

Code No. 071-06931 (10 μ g)
077-06933 (1 mg)

Granulocyte Colony-Stimulating Factor, Mouse, recombinant

顆粒球コロニー刺激因子, マウス, 組換え体
(rmG-CSF)

Source : *E. coli* expressed mouse granulocyte colony-stimulating factor

Formulation : Lyophilized from a 0.2 μ m filtered solution containing 10 mmol/L Sodium Citrate, pH 3.0.

Molecular weight : 19,100

Biological activity : ED₅₀ is below 0.05 ng/mL determined by the NFS-60 cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other aqueous buffers.

Endotoxin : Below 0.1 EU/ μ g

【Expiration date】

Indicated on the label.

【Storage】

Store at -20°C.

【Package】

Code No.	Packaging
071-06931	10 μ g
077-06933	1 mg

【Caution】

Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors

