

FUJIFILM

Wako

Code No. 077-06891 (20 μ g)
073-06893 (1 mg)

**Granulocyte-Macrophage Colony-Stimulating
Factor, Mouse, recombinant, Animal-derived-free**
顆粒球 - マクロファージコロニー刺激因子,
マウス, 組換え体, 動物由来物フリー
(rmGM-CSF, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed mouse granulocyte-macrophage colony-stimulating factor

Formulation : Lyophilized from a 0.2 μ m filtered solution containing 10 mmol/L Acetic Acid.

Molecular weight : 14,300 (calculated)

Biological activity : The ED₅₀ is less than 0.05 ng/mL determined by the FDCP-1 cell proliferation.

Reconstitution : Reconstitute in sterile water to a concentration of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other aqueous buffers.

Endotoxin : Less than 0.1 EU/ μ g

[Expiration date] Indicated on the label.

[Storage] Store at -20°C .

[Package]

Code No.	Packaging
077-06891	20 μ g
073-06893	1 mg

[Caution] Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://flwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA1