

FUJIFILM

Wako

Code No. 089-10701 (50 μ g)
085-10703 (1 mg)

**Heregulin- β -1, Human, recombinant,
Animal-derived-free**
**ヘレグリン- β -1, ヒト, 組換え体,
動物由来物フリー**

Produced using all non-animal reagents.

Source : *E. coli* expressed human heregulin- β -1

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 0.1% Trifluoroacetic Acid.

Molecular weight : 7,600 (calculated)

Biological activity : ED₅₀ is less than 0.3 ng/mL determined
by the proliferation of MCF-7 cells.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Less than 0.1 EU/ μ g.

[Expiration date] Indicated on the label.

[Storage] Store at -20°C.

[Package]

Code No.	Packaging
089-10701	50 μ g
085-10703	1 mg

[Caution] Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311-100
<http://www.wako-chemicals.de>

2309KA1