

FUJIFILM

Wako

Code No. 094-07621 (10 μ g)
090-07623 (1 mg)

Interleukin-1 β , Mouse, recombinant
インターロイキン-1 β , マウス, 組換え体
(rmlL-1 β)

Source : *E. coli* expressed mouse interleukin-1 β

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 10 mmol/L Sodium Phosphate, pH 7.5.

Molecular weight : 17,500

Biological activity : ED₅₀ is below 0.015 ng/mL determined by
the D10S cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 1 EU/ μ g

【Expiration date】 Indicated on the label.

【Storage】 Store at -20°C.

【Package】

Code No.	Packaging
094-07621	10 μ g
090-07623	1 mg

【Caution】 Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors

