

FUJIFILM

Wako

Code No. 090-07841 (100 μ g)

**Interleukin-1 Receptor Antagonist, Human,
recombinant, Animal-derived-free**
**インターロイキン-1 レセプターアンタゴニスト,
ヒト, 組換え体, 動物由来物フリー**
(rhIL-1RA, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human interleukin-1 receptor antagonist

Formulation : Lyophilized from a 0.2 μ m filtered solution containing 10 mmol/L Sodium Phosphate, pH 7.5.

Molecular weight : 17,300

Biological activity : ED₅₀ is below 50 ng/mL determined by the inhibition of human IL-1 α induced proliferation of D10.G4.1 cells assay.

Reconstitution : Reconstitute in sterile water to a concentration of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other aqueous buffers.

Endotoxin : Below 1 EU/ μ g

[Expiration date]

Indicated on the label.

[Storage]

Store at -20°C.

[Package]

100 μ g

[Caution]

Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors



2512KA1