

FUJIFILM

Wako

Code No. 096-07701 (10  $\mu$ g)  
092-07703 (1 mg)

**Interleukin-10, Human, recombinant,  
Animal-derived-free**  
**インターロイキン-10, ヒト, 組換え体,  
動物由来物フリー**  
(rhIL-10, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human interleukin-10

Formulation : Lyophilized from a 0.2  $\mu$ m filtered solution  
containing 10 mmol/L Sodium Phosphate, pH 7.5.

Molecular weight : 18,800

Biological activity : ED<sub>50</sub> is below 5 ng/mL determined by the  
MC/9 cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration  
of 0.1 mg/mL. *Do not vortex.*  
This solution can then be diluted into other  
aqueous buffers.

Endotoxin : Below 1 EU/ $\mu$ g

**[Expiration date]**

Indicated on the label.

**[Storage]**

Store at -20°C.

**[Package]**

Code No.	Packaging
096-07701	10 $\mu$ g
092-07703	1 mg

**[Caution]**

Centrifuge the vial prior to opening.  
Avoid repeated freeze-thaw cycle.

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://ffwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-311-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

2404KA1