

FUJIFILM

Wako

Code No. 097-07611 (10 μ g)
093-07613 (1 mg)

**Interleukin-3, Mouse, recombinant,
Animal-derived-free**
**インターロイキン -3, マウス, 組換え体,
動物由来物フリー**
(rmIL-3, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed mouse interleukin-3

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 0.1% Trifluoroacetic Acid.

Molecular weight : 15,200

Biological activity : ED₅₀ is below 0.25 ng/mL determined by
the NFS-60 cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 0.1 EU/ μ g

【Expiration date】

Indicated on the label.

【Storage】

Store at -20°C.

【Package】

Code No.	Packaging
097-07611	10 μ g
093-07613	1 mg

【Caution】

Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.fujifilm.co.jp>

Group Companies



Distributors

