

FUJIFILM

Wako

Code No. 094-07741 (10 μ g)
090-07743 (1 mg)

**Interleukin-1 α , Human, recombinant,
Animal-derived-free**
インターロイキン-1 α , ヒト, 組換え体, 動物由来物フリー
(rhIL-1 α , Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human interleukin-1 α

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 10 mmol/L Sodium Phosphate, pH 7.5.

Molecular weight : 18,000

Biological activity : ED₅₀ is below 0.05 ng/mL determined by
the D10.G4.1 cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 1 EU/ μ g

【Expiration date】 Indicated on the label.

【Storage】 Store at -20°C.

【Package】

Code No.	Packaging
094-07741	10 μ g
090-07743	1 mg

【Caution】 Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors



2502KA1