

FUJIFILM

Wako

Code No. 096-07681 (50 μ g)
092-07683 (1 mg)

**Insulin-like Growth Factor-I, Mouse,
recombinant**
インスリン様増殖因子-I, マウス, 組換え体
(rmlGF-I)

Source : *E. coli* expressed mouse insulin-like growth factor-I

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 0.1% Trifluoroacetic Acid.

Molecular weight : 7,700

Biological activity : ED₅₀ is below 20 ng/mL determined by the
FDC-P1 cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 0.1 EU/ μ g

【Expiration date】

Indicated on the label.

【Storage】

Store at -20°C.

【Package】

Code No.	Packaging
096-07681	50 μ g
092-07683	1 mg

【Caution】

Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2402KA1