

FUJIFILM

Wako

Code No. 116-01151 (10 µg)

KGF Solution, Human, recombinant (Keratinocyte Growth Factor)

【Source】

Corynebacterium glutamicum expressed human keratinocyte growth factor

【Formulation】

Phosphate buffer, pH 7.0 (0.2 µm filtered)

【Protein Concentration】

Indicated on the label

【Purity】

≥ 95% (by SDS-PAGE)

【Biological Activity】

The EC₅₀ is 0.1 ~ 30 ng/mL

Measurement of ERK 1/2 phosphorylation assay using MCF-7 cells

【Endotoxin】

< 0.1 EU/µg

【Storage】

Store at -80°C

Avoid repeating freeze-thaw.

【Package】

Code No.	Package
116-01151	10 µg

【Caution】

Centrifuge the vial prior to opening.

【Description】

This product is manufactured by Ajinomoto Co., Inc.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>