

FUJIFILM

Wako

Code No. 120-05751 (1 mg)

Laminin Solution, from Mouse EHS Tumor
ラミニン溶液, マウスEHS肉腫由来

Laminin is an extracellular matrix protein and a major structural component of basement membranes. Laminin is used for the coating of culture dishes. Laminin mainly promotes the attachment of epithelial cells, muscle cells, hepatocytes and so on.

This product is for laboratory use only ; use in any such application is the responsibility of the user.

Source : Mouse Engelbreth-Holm-Swarm Tumor

Formulation : Sterile-filtered 50 mmol/L Tris-HCl (pH7.4), 150 mmol/L NaCl

Protein Concentration: Indicated on the label

Storage : Store at -80°C
Avoid repeated freeze and thaw cycles. When stored at -20°C , this product is stable for a year.

Package : 1 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggenstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA2