

FUJIFILM

Wako

Code No. 121-07001 (1 mg)

**Leptin, Mouse, recombinant,  
Animal-derived-free**  
**レプチン, マウス, 組換え体, 動物由来物フリー**  
(rmOB, Animal-derived-free)

---

Produced using all non-animal reagents.

Source : *E. coli* expressed mouse leptin

Formulation : Lyophilized from a 0.2  $\mu$ m filtered solution  
containing 10 mmol/L Sodium Phosphate, pH 7.5.

Molecular weight : 16,100

Reconstitution : Reconstitute in sterile water to a concentration  
of 0.1 mg/mL. *Do not vortex.*

This solution can then be diluted into other  
aqueous buffers.

Endotoxin : Below 1 EU/ $\mu$ g

**【Expiration date】** Indicated on the label.

**【Storage】** Store at -20°C.

**【Package】** 1 mg

**【Caution】** Centrifuge the vial prior to opening.  
Avoid repeated freeze-thaw cycle.

---

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors

