

**FUJIFILM**

**Wako**

< for Research Use Only >

Code No. 127-06361 (1 mg)

**LSL-N Lectin, recombinant, Solution**  
**LSL-N レクチン, 組換え体, 溶液**

LSL protein is composed of the N-terminal domain known as the hemolytic factor and polysaccharide binding factor and the C-terminal domain known as the pore-forming factor. LSL-N lectin is expressed as the N-terminal domain.

**Source :** *E. coli* expressed *Laetiporus sulphureus* lectin (LSL-N).

**Appearance :** PBS (pH 7.4).

**Protein concentration :** Indicated on the label.

**Glycan specificity :** Poly-LacNAc

**Manufacture :** Affinity purification by lactose agarose.

**[Expiration date]** Indicated on the label.

**[Storage]** Store at  $-20^{\circ}\text{C}$ . After thawing, keep at  $2-10^{\circ}\text{C}$ .  
Avoid the freeze-thaw cycles.

**[Package]**

| Code No.  | Packaging |
|-----------|-----------|
| 127-06361 | 1 mg      |

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-3111-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

2305KA2