

FUJIFILM

Wako

< for Research Use Only >

Code No. 136-18221 (1 mg)
132-18223 (1 mg × 5)

Malectin Lectin, recombinant, Solution
マレクチンレクチン, 組換え体, 溶液

Malectin is an endoplasmic reticulum resident lectin involving in glycoprotein quality control and known that it is highly conserved in animals. It is known that specially bind to Glc₂-high-mannose *N*-glycan.

Source : *E. coli* expressed human malectin.

Appearance : PBS (pH 7.4).

Protein concentration : Indicated on the label.

Carbohydrate specificity : Glc₂-*N*-glycan

Manufacture : Affinity purification by maltose agarose.

[**Expiration date**] Indicated on the label.

[**Storage**] Store at -20 °C. After thawing, keep at 2-10 °C.
Avoid the freeze-thaw cycles.

[**Package**]

Code No.	Packaging
136-18221	1 mg
132-18223	1 mg × 5 vials

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2