

FUJIFILM

Wako

< for Research Use Only >

Code No. 137-18011 (1 mg)

MOA Lectin, recombinant, Solution
MOAレクチン, 組換え体, 溶液

MOA lectin is isolated from the fruiting bodies of the mushroom *Marasmius oreades*. It is known that specifically bind to non-reducing terminal Gal α 1-3Gal-carbohydrates.

Source : *E. coli* expressed *Marasmius oreades* lectin (MOA).

Appearance : 0.1 $\times$ PBS (pH 7.4).

Protein concentration : Indicated on the label.

Glycan specificity : Gal α 1-3Gal
Gal α 1-3Gal β 1-4GlcNAc
Gal α 1-3 [Fuc α 1-2] Gal β 1-4GlcNAc
[Blood Group B]

Manufacture : Affinity purification by melibiose agarose.

[Expiration date] Indicated on the label.

[Storage] Store at -20°C . After thawing, keep at $2-10^{\circ}\text{C}$.
Avoid the freeze-thaw cycles.

[Package]

Code No.	Packaging
137-18011	1 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2