

FUJIFILM

Wako

Code No. 139-19551 (10 μ g)
135-19553 (1 mg)

**Myostatin, Human, recombinant,
Animal-derived-free**
**ミオスタチン, ヒト, 組換え体,
動物由来物フリー**
(rhMSTN/GDF-8, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human myostatin

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 0.1% Trifluoroacetic Acid.

Molecular weight : 24,800 (dimer)

Biological activity : ED₅₀ is below 50 ng/mL determined by the
MPC-11 cell cytotoxicity assay.

Reconstitution : Reconstitute in sterile 20 mmol/L HCl at 0.1
mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 0.1 EU/ μ g

【Expiration date】

Indicated on the label.

【Storage】

Store at -20°C.

【Package】

Code No.	Packaging
139-19551	10 μ g
135-19553	1 mg

【Caution】

Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2402KA1