

FUJIFILM

Wako

Code No. 147-09521 (10 mg)
143-09523 (10 mg × 5)

Normal Rat IgG, Whole Molecule, Purified
正常ラットIgG, 全分子, 精製品

This product is purified from normal rat serum through chromatography. It can be used as control, standard, blocking agent, or coating protein in a variety of assays, including ELISA, immunoblotting, immunoprecipitation, immunodiffusion and immunoelectrophoresis.

Source Normal Rat Serum

Form Supplied in 10 mmol/L Sodium Phosphate, 0.15 mol/L Sodium Chloride, pH 7.2 with 0.05w/v% Sodium Azide.

Protein Concentration Indicated on the label.

Purification Protein G affinity chromatography

Storage 2-10 °C

Package

Code No.	Packaging
147-09521	10 mg
143-09523	10 mg × 5

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggenstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA2