

**FUJIFILM**

**Wako**

Code No. 141-09541 (10 mg)  
147-09543 (10 mg × 5)

**Normal Bovine IgG, Whole Molecule, Purified**  
**正常ウシIgG, 全分子, 精製品**

This product is purified from normal bovine serum through chromatography. It can be used as control, standard, blocking agent, or coating protein in a variety of assays, including ELISA, immunoblotting, immunoprecipitation, immunodiffusion and immunoelectrophoresis.

**Source** Normal Bovine Serum

**Form** Supplied in 10 mmol/L Sodium Phosphate, 0.15 mol/L Sodium Chloride, pH 7.2 with 0.05w/v% Sodium Azide.

**Protein Concentration** Indicated on the label.

**Purification** Protein A affinity chromatography

**Storage** 2-10 °C

**Package**

| Code No.  | Packaging |
|-----------|-----------|
| 141-09541 | 10 mg     |
| 147-09543 | 10 mg × 5 |

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://ffwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggenstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-3111-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

2306KA2