

FUJIFILM**Wako**

Code No. 148-09551 (10 mg)

Normal Rabbit IgG, Whole Molecule, Purified
正常ウサギIgG, 全分子, 精製品

This product is purified from normal rabbit serum through chromatography. It can be used as control, standard, blocking agent, or coating protein in a variety of assays, including ELISA, immunoblotting, immunoprecipitation, immunodiffusion and immunoelectrophoresis.

Source Normal Rabbit Serum**Form** Supplied in 10 mmol/L Sodium Phosphate, 0.15 mol/L Sodium Chloride, pH 7.2 with 0.05w/v% Sodium Azide.**Protein Concentration** Indicated on the label.**Purification** Protein A affinity chromatography**Storage** 2-10 °C**Package**

Code No.	Packaging
148-09551	10 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggenstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA2