

FUJIFILM**Wako**

Code No. 143-09501 (10 mg)
149-09503 (10 mg × 5)

Normal Human IgG, Whole Molecule, Purified
正常ヒトIgG, 全分子, 精製品

This product is purified from normal human serum/plasma through chromatography. It can be used as control, standard, blocking agent, or coating protein in a variety of assays, including ELISA, immunoblotting, immunoprecipitation, immunodiffusion and immunoelectrophoresis.

Source	Normal Human Serum/Plasma
Form	Supplied in 10 mmol/L Sodium Phosphate, 0.15 mol/L Sodium Chloride, pH 7.2 with 0.05w/v% Sodium Azide.
Protein Concentration	Indicated on the label.
Purification	Protein A affinity chromatography
Disease Testing	Source serum was found negative for antibodies to HIV 1/2 and HCV, and for HBsAg, HIV-1 RNA, HCV RNA and Syphilis.
BIOHAZARD	Handle as if capable of transmitting infectious agents.
Storage	2-10℃

Package

Code No.	Packaging
143-09501	10 mg
149-09503	10 mg × 5

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA2