

FUJIFILM

Wako

Code No. 160-19501 (1 vial/1 mL)

**Protease Inhibitor Mixture-DMSO Solution
for Mammalian Cell and Tissue Extracts**

プロテアーゼ阻害剤ミックス, DMSO 溶液
(動物細胞・組織抽出液用)

Appearance : DMSO solution (frozen)

Contents :

Inhibitor	M. W.	Concentration in the vial	Target protease
AEBSF Hydrochloride	239.70	100 mmol/L	Serine Proteases
Aprotinin, Bovine Lung	6511.45	80 μ mol/L	Serine Proteases (Broad Spectrum)
Bestatin	308.38	5 mmol/L	Aminopeptidase B Leucine Aminopeptidase
E-64, Protease Inhibitor	357.41	1.5 mmol/L	Cysteine Proteases
Leupeptin Hemisulfate	475.59	2 mmol/L	Cysteine Proteases Trypsin-like Protease
Pepstatin A	685.89	1 mmol/L	Aspartic Proteases

Storage : Store at -20°C

Package : 1 vial (1mL)

Precaution : Handle material with protective glasses,
protective gloves and protective mask.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggenstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-3111100
<http://www.wako-chemicals.de>

2305KA2