

FUJIFILM**Wako**

Code No. 160-24371 (1 mL)

**Phosphatase Inhibitor Cocktail
Solution II (× 100)****ホスファターゼ阻害剤カクテル溶液Ⅱ (× 100)**

Phosphatase Inhibitor Cocktail Solution II contains five phosphatase inhibitors that will inhibit acid and alkaline phosphatases as well as protein tyrosine phosphatases. This product is used for cell lysates and tissue extracts.

This product is for laboratory use only ; use in any such application is the responsibility of the user.

Formulation : Aqueous solution**Components** :

Phosphatase Inhibitor	M.W.(anhyd.)	Target Phosphatase
Sodium Fluoride	41.99	Acid Phosphatases
Sodium Orthovanadate(V)	183.91	Tyrosine Phosphatases, Alkaline Phosphatases
Sodium Tartrate	194.05	Acid Phosphatases
Imidazole	68.08	Alkaline Phosphatases
Sodium Molybdate(IV)	205.94	Acid Phosphatases

Storage : Store at 2~10℃**Package** : 1 mL**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2304KA2