

FUJIFILM

Wako

Code No. 166-29231 (50 μ g)
162-29233 (1 mg)

Prolactin, Human, recombinant プロラクチン, ヒト, 組換え体

Source : *E. coli* expressed human prolactin

Formulation : Lyophilized from a 0.2 μ m filtered solution containing 10 mmol/L Sodium Phosphate, pH 7.5.

Molecular weight : 23,000 (calculated)

Biological activity : ED₅₀ is less than 1 ng/mL determined by the NB2-11 cell proliferation.

Reconstitution : Reconstitute in sterile water to a concentration of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other aqueous buffers.

Endotoxin : Less than 1 EU/ μ g

[Expiration date]

Indicated on the label.

[Storage]

Store at -20°C.

[Package]

Code No.	Packaging
166-29231	50 μ g
162-29233	1 mg

[Caution]

Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA1